



Please Read Before Use

The Vicair® Liberty Air cushion is an adjustable cushion. The cushion itself has a set of instructions printed on the bottom panel that clearly indicate it's proper setup. After the cushion has been properly adjusted for the user, please be sure to close and recess the blue valve into the cushion.

An optional cover has been supplied in the separate bag and can be installed if desired. To install the cover properly be sure to align the blue valve of the cushion to the hole located in the cover for future access. A set of instructions is included in the bag with the cover.

Note that the zipper side goes down and the valve portion is always located to the right rear of the patient.

Toll Free 800.564.9248 www.comfortcompany.com



IS-PCLIBERTYPREUSE
REV0413



Please Read Before Use

The Vicair® Liberty Air cushion is an adjustable cushion. The cushion itself has a set of instructions printed on the bottom panel that clearly indicate it's proper setup. After the cushion has been properly adjusted for the user, please be sure to close and recess the blue valve into the cushion.

An optional cover has been supplied in the separate bag and can be installed if desired. To install the cover properly be sure to align the blue valve of the cushion to the hole located in the cover for future access. A set of instructions is included in the bag with the cover.

Note that the zipper side goes down and the valve portion is always located to the right rear of the patient.

Toll Free 800.564.9248 www.comfortcompany.com



IS-PCLIBERTYPREUSE
REV0413