UNIVERSAL HALO SAFETY RING

PRODUCT GUIDE

| PART NUMBER | DESCRIPTION | RETAIL |
|-------------|----------------------------|----------|
| 77301 | Halo Assembly with 1 Ring | \$387.50 |
| 77302 | Halo Assembly with 2 Rings | \$600.00 |

^{* 77101} has been discontinued and replaced by 77301

IMPROVEMENTS & UPDATES

The new Halo assembly fits on more beds including all of the Joerns Easy Care Line and Invacare CS Line.

New assembly replaces 77113 and 77203 that were unable to articulate with the bed frame.

Both the 77301 and 77302 come with 2 ring mount brackets and able to mount on either side of the bed for single assemblies.


New assembly hardware includes center bar with mounting holes, screws, and 2 Allen wrenches.

Initial adjustments are based on the bed width and can be done prior to installing on frame using 1 Allen wrench.

Final adjustments are done at each side of the bed frame using 1 Allen wrench.

Also includes cone tip screws that vertically lock into the bed frame to help stabilize and secure the Halo assembly.

The new assembly has improved strength and durability.



UNIVERSAL HALO FIT LIST

| HOSPITAL BED | PART NUMBER | |
|----------------------------|-----------------------------------|--|
| Panacea 1000 | 77301/77302 | |
| Panacea 1500 | 77401/77402 | |
| Panacea 6300 Bariatric | 77301/77302 | |
| Joerns Care 100 35" W | 77301/77302 | |
| Joerns UCXT 35" W | 77301/77302 | |
| Joerns EasyCare Line 35" W | 77301/77302 | |
| Joerns EasyCare Line 39" W | 77301/77302 | |
| Joerns EasyCare Line 42" W | 77301/77302 | |
| Invacare CS Line 36" W | 77301/77302 | |
| Invacare CS Line 39" W | 77301/77302 | |
| Invacare CS Line 42" W | 77301/77302 | |
| Invacare CS9 FX600 | Unable to Attach to Halo Products | |
| Noa Elite Riser 35" W | Unable to Attach to Halo Products | |

| HOSPITAL BED | PART NUMBER |
|--|-----------------------------------|
| Joerns EC2000 (discont.) | 77301/77302 |
| Invacare SC900 | 77301/77302 |
| Invacare 820 | 77301/77302 |
| Invacare Carol (Arro, Echo, Spirit) | 77103* |
| Invacare CS C59FX600 36"-42" W | 77301/77302 |
| Simmons (All) | 77301/77302 |
| Medline 1232 | Unable to Attach to Halo Products |
| Zenith PBH093A35MA | Unable to Attach to Halo Products |

^{*}Universal Halo attaches to additional bed models not listed above. *77103 = Carol Mounting Hardware.

^{* 77102} has been discontinued and replaced by 77302