

# QUICK RELEASE AMP

**E1020**  
RESIDUAL LIMB SUPPORT SYSTEM



## PRODUCT FEATURES

- Amp support may be used with any wheelchair with a solid seat pan.
- Installs onto either left or right side of solid seat pan.
- Toggle lever allows positioning at multiple angles.
- Amputee pad has adduction to help center residual limb.
- Foam construction has soft entry layer for comfort and firmer base layer for support.
- TB3 parts designed for use on Quantum Rehab Tru-Balance® 3 chairs.

## QUICK RELEASE AMP HARDWARE

Fast and simple removal with the quick release assembly.



## ANGLE ADJUSTMENT

Toggle lever allows positioning at multiple angles.



### 10" QUICK-RELEASE AMP PART NUMBERS

W x D	Thickness	Extension Kit
8" x 10"	3"	AMPKIT810
9" x 10"	3"	AMPKIT910
10" x 10"	3"	AMPKIT1010
11" x 10"	3"	AMPKIT1110
12" x 10"	3"	AMPKIT1210

### 10" TB3 QUICK-RELEASE AMP PART NUMBERS

W x D	Thickness	Extension Kit
8" x 10"	3"	AMPTQR810
9" x 10"	3"	AMPTQR910
10" x 10"	3"	AMPTQR1010
11" x 10"	3"	AMPTQR1110
12" x 10"	3"	AMPTQR1210

### 14" QUICK-RELEASE AMP PART NUMBERS

W x D	Thickness	Extension Kit
8" x 14"	3"	AMPKIT814
9" x 14"	3"	AMPKIT914
10" x 14"	3"	AMPKIT1014
11" x 14"	3"	AMPKIT1114
12" x 14"	3"	AMPKIT1214

### 14" TB3 QUICK-RELEASE AMP PART NUMBERS

W x D	Thickness	Extension Kit
8" x 14"	3"	AMPTQR814
9" x 14"	3"	AMPTQR914
10" x 14"	3"	AMPTQR1014
11" x 14"	3"	AMPTQR1114
12" x 14"	3"	AMPTQR1214

Comes standard with Comfort-Tek. See Pg. 2 for details.