

SWING-AWAY AMP

E1020
RESIDUAL LIMB SUPPORT SYSTEM



PRODUCT FEATURES

- Designed to install on all popular wheelchairs within minutes.
- Provides patients a fully adjustable amputee support that will coincide with any wheelchair cushion.
- Allows users to be transferred in and out of a wheelchair within seconds.
- Installs onto either left or right side of wheelchair.
- Fits tubing sizes from 7/8" to 1".
- Amputee support pad has adduction to help center residual limb.

10" D SWING-AWAY AMP PART NUMBERS

W x D	Thickness	Extension Kit
8" x 10"	3"	AMPSA810
9" x 10"	3"	AMPSA910
10" x 10"	3"	AMPSA1010
11" x 10"	3"	AMPSA1110
12" x 10"	3"	AMPSA1210

14" D SWING-AWAY AMP PART NUMBERS

W x D	Thickness	Extension Kit
8" x 14"	3"	AMPSA814
9" x 14"	3"	AMPSA914
10" x 14"	3"	AMPSA1014
11" x 14"	3"	AMPSA1114
12" x 14"	3"	AMPSA1214

Comes standard with Comfort-Tek. See Pg. 2 for details.

SWING-AWAY AMP HARDWARE

U.S. Patent No. 8,231,090

Hardware allows for fast and simple swing away and removal.

Swing away feature allows patients to be easily transferred in and out of a wheelchair.

Pad can be easily removed from and reinstalled onto wheelchair.



ANGLE ADJUSTMENT

Toggle lever allows positioning at multiple angles.

Up to 4" of height adjustment to accommodate varying seat cushion thicknesses.

Pad can be adjusted 3" in depth to accommodate wheelchair depth conversion kits.

